



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Anesthesia machine with LED Screen

Product Code :
SAMCOMEDICAL3001A054

Description :

Anesthesia machine with LED Screen

Technical Specification :

Anesthesia machine with LED Screen With LED screen. Mounted on 4 castors. 2 front castors with brakes. With 2 big drawers . Handle for easy maneuvering Epoxy painted high quality mild steel frame. Equipped with Anesthesia Vaporizer. Agent type: Isoflurane, Enflurane, Halothane. Classical CO2 circle absorber with 1 canister With necessary monitoring and alarm function, to ensure the operational safety



Samco Medical

Website: www.samcomedical.com,