



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Anesthesia machine with 5.5 SNT screen

Product Code :

SAMCOMEDICAL3001A055

Description :

Anesthesia machine with 5.5 SNT screen

Technical Specification :

Anesthesia machine with 5.5 SNT screen With 5.5" SNT screen. Mounted on 4 castors. 2 front castors with brakes. With 2 big drawers. Handle for easy maneuvering. Epoxy painted high quality mild steel frame. Equipped with High precise Anesthesia Vaporizer (with temperature, flow and pressure compensation function). The agent output is continuous, stable & precise. The agent type: Isoflurane, Enflurane, Halothane, Sevoflurane Anesthesia ventilator with IPPV and SIMV mode. The ventilator can display P-t and V-t wave-form at the same time. Necessary monitoring and alarm function to ensure operational safety. Classical CO2 circle absorber with 2 canisters for low resistance and minimum gas leakage. Optional : O2, N2O Yoke device. When main gas supply fails, the backup gas supply starts



Samco Medical

Website: www.samcomedical.com,