



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Anesthesia Machines with CO2, Multi Parameter,  
Anesthesia Depth Monitor

**Product Code :**

SAMCOMEDICAL3001A056

**Description :**

Anesthesia Machines with CO2, Multi Parameter, Anesthesia Depth Monitor

**Technical Specification :**

Anesthesia Machines with CO2, Multi Parameter, Anesthesia Depth Monitor With 5.5 SNT screen. Mounted on 4 castors. 2 front castors with brakes. With 1 safety lock and 2 big drawers Handle for easy maneuvering. Aluminium alloy rail frame. With durable and non-corrosive stainless steel platform. Equipped with high precise Anesthesia Vaporizer (with temperature, flow and pressure compensation function). The agent output is continuous, stable & precise. The agent type: Isoflurane, Enflurane. With independent Anesthesia Ventilator, the ventilator frame is moulded and height can be adjusted. The ventilator can display P-t and V-t wave-form at the same time. Classical CO2 circle absorber with 2 canisters for low resistance and minimum gas leakage. Fresh gas outlet is in front of frame for the user to easily connect Jackson-Rees circuit. With upgraded functions, eg, CO2 Monitor, Multi-parameter monitor, Anesthesia depth monitor, etc. Optional : O2, N2O Yoke device. When main gas supply is fails, the backup gas supply starts.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),