



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Anesthesia Machine with 10.4 TFT screen

Product Code :

SAMCOMEDICAL3001A057

Description :

Anesthesia Machine with 10.4 TFT screen

Technical Specification :

Anesthesia Machine with 10.4 TFT screen With 10.4 TFT color screen. Mounted on 4 castors. 2 front castors with brakes. With 2 big drawers. Handle for easy maneuvering. Epoxy painted high quality mild steel frame. Equipped with Two high precise Anesthesia Vaporizer (with temperature, flow and pressure compensation function). The agent output is continuous, stable & precise. The agent type: Isoflurane, Enflurane, Halothane, Sevoflurane. Classical CO2 circle absorber with 2 canisters for low resistance and minimum gas leakage. Optional: O2, N2O Yoke device. When main gas supply is fails, the Backup gas supply starts.



Samco Medical

Website: www.samcomedical.com,