



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Anesthesia machine integrated for Pediatric and Adult

**Product Code :**

SAMCOMEDICAL3001A058

**Description :**

Anesthesia machine integrated for Pediatric and Adult

**Technical Specification :**

Anesthesia machine integrated for Pediatric and Adult With 10.4 TFT color screen Mounted on 4 castors. 2 front castors with brakes. With 1 safety lock and 3 big drawers. Handle for easy maneuvering. Aluminium alloy rail frame. With durable and non-corrosive stainless steel platform. Equipped with high precise Anesthesia Vaporizer (with temperature, flow and pressure compensation function). The agent output is continuous, stable & precise The agent type: Isoflurane, Enflurane, Halothane, Sevoflurane. With Volume & pressure ventilation modes Classical CO2 circle absorber with 2 canisters for low resistance and minimum gas leakage Fresh gas outlet is in front of frame for the user to easily connect Jackson-Rees circuit. With upgraded with functions, eg, CO2 Monitor, Multi-parameter monitor, Anesthesia depth monitor, etc Working interfaces for adult and pediatric are separated; the doctor can choose interface according to patient type Optional: O2, N22 Yoke device. When main gas supply is fails, The backup gas supply starts.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),