



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Ceiling Surgical Operating Light Small Single

**Product Code :**  
SAMCOMEDICAL200001A012



**Description :**

Ceiling Surgical Operating Light Small Single has daylight characters and can be adjusted as per the operation requirements. Lights arms are made of Stainless Steel Fitting. The range of the color temp - 4,200 K. Detachable handle for simplistic sterilization. Parabolic reflectors & color correcting filters. Dome Head: 550 mm Halogen Bulbs: 24V x 70 W x 4 nos. Light intensity at 1 mt: 50,000 Lux Height Adjustment: 914 mm Parabolic Reflector: 180 mm Color Temperature: 4200 K Control Unit: Transformer Ceiling Surgical Operating Light Small Single, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),