



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Rack Steam Sterilizer, Double Rack

Product Code :

SAMCOMEDICAL5006A017

Description :

Rack Steam Sterilizer, Double Rack

Technical Specification :

Rack Steam Sterilizer, Double Rack Steam sterilizer, pressure type, double rack, used to sterilize reusable: syringes and needles. Aluminium vessel with high pressure seal suitable for carrying out sterilization using superheated steam. Steam sterilizer designed to operate on external heat source (non-electric) at $+121^{\circ}\pm\frac{1}{2}^{\circ}\text{C}$ over a 20 minute period. Components: Aluminium vessel with handles. Aluminium cover with handles and rubber sealing gasket. Pressure valve. Safety plug. 2 aluminium racks with 2 lid, for 64 syringes (44 syringes 2 ml, 16 syringes 1 ml, 4 syringes 5 ml) and 80 needles per rack. Sterilization capacity: 9 litres. Supplied with: 1 x sterilization indicator sticker (pack of 100), 2 x spare safety plugs, 3 x spare sealing gaskets, 1 x plastic bowl for washing syringes, 1 x carrying bag, 1 x user manual in English.



Samco Medical

Website: www.samcomedical.com,