



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Autoclaves Pressure Steam Sterilizers (Aluminum) -  
Economy Mirror / Polished Finish

**Product Code :**

SAMCOMEDICAL5006A019

**Description :**

Autoclaves Pressure Steam Sterilizers (Aluminum) - Economy Mirror / Polished Finish

**Technical Specification :**

Autoclaves Pressure Steam Sterilizers (Aluminum) - Economy Mirror / Polished Finish Construction: These autoclaves are made of 10 gauge aluminum. The drum is deep drawn & the ring on top is riveted. The lid is made by casting. Locking of lid by 6 nuts with rubber gasket seal in between (not metal to metal seal). Two safety valves, pressure gauge, flexible metal rubbing to exhaust air. Includes outer stand, inner aluminum container, inner tripod stand.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),