



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Autoclaves Pressure Steam Sterilizers (Stainless Steel) - Seamless

Product Code :

SAMCOMEDICAL5006A021

Description :

Autoclaves Pressure Steam Sterilizers (Stainless Steel) - Seamless

Technical Specification :

Autoclaves Pressure Steam Sterilizers (Stainless Steel) - Seamless Construction: These portable vertical autoclaves are made out of thick high quality (14 gauge) stainless steel sheet, deep drawn to cylindrical shape, having no joints and welding. Seamless construction allows no bacteria residue and free from dirt accumulation. Fitted with pressure gauge, safety valve (weight type), seam/pressure release cock, lid handle. Lid is similar to pressure cooker with locking over rubber gasket Features: Stainless steel seamless construction. Portable vertical pressure type, easy to use top loading autoclaves. The steam pressure inside with autoclave can develop from 1.10 kgf / sq.cm (15 psi) to 1.25 kgf / sq.cm (18 psi). Equipped with top handle pressure gauge and pressure release valve.



Samco Medical

Website: www.samcomedical.com,