



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Light Weight Disinfection Tray System

Product Code :
SAMCOMEDICAL5006A030

Description :

Light Weight Disinfection Tray System

Technical Specification :

Light Weight Disinfection Tray System This tray system is a nice sterilizing / disinfecting / soaking tray that is essential for sterilization of certain instruments & equipment which cannot be sterilized under steam & pressure in Autoclaves. These can be also used as a back up to your other sterilizing system. The tray includes Outer Tray with lid & a Lift-out drain basket with handles to facilitate lifting of instruments in and out of the soaking solution. Lid protects against spillage and minimizes vaporization Dimension: Outer tray measures 845mm L x 225mmWx 150mm H Volume Capacity: 10 liters appx, [NON AUTOCLAVABLE]



Samco Medical

Website: www.samcomedical.com,