



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Double Wall Fully Automatic Autoclave

Product Code :
SAMCOMEDICAL5006A031

Description :

Double Wall Fully Automatic Autoclave

Technical Specification :

Double Wall Fully Automatic Autoclave Top end model for sterilization of all critical application with exact reproducibility Unique double wall design Heavy duty jacket Fully automatic operation $\frac{1}{2}$ P based PID controller with dual display for Temp. & Time Pre & Post Vacuum available Excellent performance in drying Reduced loading height



Samco Medical

Website: www.samcomedical.com,