



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Table Top Front Loading Autoclave Round Chamber

Product Code :

SAMCOMEDICAL5006A034

Description :

Table Top Front Loading Autoclave Round Chamber

Technical Specification :

Table Top Front Loading Autoclave Round Chamber Table top front loading autoclave Fully Automatic micro-processor based control Sterilization chamber made up of thick stainless steel sheet, tested for high pressure Digital timer for both wet and dry cycle Auto drain of water to reservoir tank & also condensation of steam Alarm after the total command is completed Cycle time 10-15 minutes



Samco Medical

Website: www.samcomedical.com,