



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Samcomedical Gas cum Electric Autoclaves

**Product Code :**

SAMCOMEDICAL5006A035

**Description :**

Samcomedical Gas cum Electric Autoclaves

**Technical Specification :**

Gas cum Electric Autoclaves Energy operated LPG operated autoclaveBiggest advantage where electricity is a major problem Steam can be generated by LPG or Electricity Vacuum drying available



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),