



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Autoclave Steam Indicator Tape

**Product Code :**  
SAMCOMEDICAL5006A037

**Description :**

Autoclave Steam Indicator Tape

**Technical Specification :**

Autoclave Steam Indicator Tape Light Yellow indicator tapes has stripes that changes colour from yellow to bark brown/black when passes through steam sterilization Cycle (i.e. 20 minutes of exposure to a temperature of 121°C or 5 minutes of exposure to a temperature of 134°C)



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),