



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Electric - S.S. 304 Grade

Product Code :
SAMCOMEDICAL5006A038

Description :

Electric - S.S. 304 Grade

Technical Specification :

Electric - S.S. 304 Grade Instrument Sterilizers, electric, with hinged lid fitted with immersion rod and water outlet tap, seamless stainless steel construction. Fitted with automatic cutoff safety pin device, thermostat & tray lifting arrangement



Samco Medical

Website: www.samcomedical.com,