



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Portable Autoclave - Semi Automatic

Product Code :
SAMCOMEDICAL5006A039

Description :

Portable Autoclave - Semi Automatic

Technical Specification :

Portable Autoclave - Semi Automatic SS Seamless die pressed deep drawn chamber ?p based PID controller with dual display for Temp. & Time Fully SS construction Dual safety valves Microprocessor Digital temperature controller cum timer Over heat auto shut down at 135 \pm 1/2C Dual Display Unique shape heater of outer SS tubing Food grade silicon rubber gasket for better durability Specially designed heater for optimum chamber usage Calibrated temperature-cum-pressure gauge Colour coded monitoring gauge SS Safety valve Vacuum breaker to remove air pockets in the chamber



Samco Medical

Website: www.samcomedical.com,