



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hospital Platelet Incubator with Agitator

Product Code :
SAMCOMEDICAL230001A002



Description :

Hospital Platelet Incubator with Agitator - Platelet Incubator with Agitator: Microprocessor based PID temperature controller for controlling accuracy of $\pm 0.1^\circ\text{C}$. Platelet incubator with agitator provide smooth platelet agitation of 70 ± 2 RPM to increase platelet viability. Forced air cooling unit in platelet incubator with agitator maintains uniform temperature. Platelet incubator with agitator comes with microprocessor based temperature indicator to control resolution to 0.1°C . Platelet incubator with agitator also has 7 days temperature chart recorder. Hospital Platelet Incubator with Agitator, Blood Bank Equipments Exporters, Blood Bank Equipments Manufacturers



Samco Medical

Website: www.samcomedical.com,