



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Hot Air Sterilizer upto 250C (Laboratory Electric Oven Universal Type)

Product Code :

SAMCOMEDICAL5006A045

Description :

Hot Air Sterilizer upto 250C (Laboratory Electric Oven Universal Type)

Technical Specification :

Hot Air Sterilizer - upto 250 \pm 1/2 $^{\circ}$ C (Laboratory Electric Oven Universal Type) Temperature Thermostatically controlled by hydraulic thermostat from 50 \pm 1/2 $^{\circ}$ C to 250 \pm 1/2 $^{\circ}$ C with a sensitivity of \pm 1 \pm 1/2 $^{\circ}$ C in Bare Oven. Beaded heating elements are placed in ribs at the bottom and sides. Built in L-shape thermometer. Double walled, inside made of Aluminium or Stainless Steel & outside mild steel sheet painted with attractive stove enamel with adjustable wire mesh nickel plated trays, on/off Rocker switch with Indicator & Thermostat Indicator, to work on 220/230 Volts A.C



Samco Medical

Website: www.samcomedical.com,