



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Hot Air Sterilizer

**Product Code :**  
SAMCOMEDICAL5006A046

**Description :**

Hot Air Sterilizer (Laboratory Electric Oven Universal Type), Triple Wall with Air Circulating System upto 250°C

**Technical Specification :**

Hot Air Sterilizer (Laboratory Electric Oven Universal Type), Triple Wall with Air Circulating System upto 250½°C In this system there is one path of circulating Hot Air in between the inner chamber and insulation. Forced air, moved by the motorized blower results in minimum temperature variation at any point in working space Temperature controlled by Digital Indicator cum Controller PID based with Safety Alarm from 50½°C to 250½°C with a Sensitivity of +/- 1½°C or better in Bare Oven. Working chamber is made of Stainless steel 304 grade and outside made of mild steel painted with good stoving enamel with adjustable wire mesh Nickel plated trays. To work on 220/230 Volts A.C. supply.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),