



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
HIGH Temperature Furnace of 1200C

**Product Code :**  
SAMCOMEDICAL5006A048

**Description :**

HIGH Temperature Furnace of 1200C

**Technical Specification :**

HIGH Temperature Furnace of 1200 $\pm$ 1/2 $\circ$ C  
Light weight with ceramic fibre wool insulation. The outer casing is made of double walled thick P.C.R.C. sheet, painted with attractive stove enamel. Heating Elements are made of KANTHAL "A-1" wire in coil Form placed in the ribs of the tiles, backed by cerawool insulation, which avoids loss of energy. A viewing cap is also provided on the door. Temperature controlled by Digital Indicator cum Controller PID based with safety Alarm an accuracy of  $\pm$  1/2  $\pm$  3 $\pm$  1/2 $\circ$ C or better. To work on 220/230 Volts or 400/440 Volts A.C. 3 phase supply. Maximum Temp. 1200 $\pm$  1/2 $\circ$ C working temp. 1150 $\pm$  1/2 $\circ$ C



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),