



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
HIGH Temperature Furnace of 1200C

Product Code :
SAMCOMEDICAL5006A048

Description :

HIGH Temperature Furnace of 1200C

Technical Specification :

HIGH Temperature Furnace of 1200 \pm 1/2C
Light weight with ceramic fibre wool insulation. The outer casing is made of double walled thick P.C.R.C. sheet, painted with attractive stove enamel. Heating Elements are made of KANTHAL "A-1" wire in coil Form placed in the ribs of the tiles, backed by cerawool insulation, which avoids loss of energy. A viewing cap is also provided on the door. Temperature controlled by Digital Indicator cum Controller PID based with safety Alarm an accuracy of \pm 1/2 3 \pm 1/2C or better. To work on 220/230 Volts or 400/440 Volts A.C. 3 phase supply. Maximum Temp. 1200 \pm 1/2C working temp. 1150 \pm 1/2C



Samco Medical

Website: www.samcomedical.com,