



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Laboratory Rectangular Heating Plate with Cast Iron Top

Product Code :

SAMCOMEDICAL5006A049

Description :

Laboratory Rectangular Heating Plate with Cast Iron Top

Technical Specification :

Laboratory Rectangular Heating Plate with Cast Iron Top This is very useful for continuous heating upto 350½C. The body is made of specially designed with perforated P.C.R.C. sheet painted with attractive stove enamel.



Samco Medical

Website: www.samcomedical.com,