



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Laboratory Round Heating Plate

**Product Code :**  
SAMCOMEDICAL5006A050

**Description :**

Laboratory Round Heating Plate

**Technical Specification :**

Laboratory Round Heating Plate The Body is made of P.C.R.C. sheet painted with attractive stove enamel. The top of plate is 20.5 cm dia. fitted with an indicating lamp. This type of heating plate is very useful for fast heating upto 500C. To work on 220/230 Volts A.C



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),