



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Laboratory Autoclaves - Semi Automatic Series

Product Code :
SAMCOMEDICAL5006A051

Description :

Laboratory Autoclaves - Semi Automatic Series

Technical Specification :

Laboratory Autoclaves - Semi Automatic Series Microprocessor based PID controller with dual display for Temp. & Time Unique pedal free design Reduced loading height No hinge or shaking lids Fully SS 304 construction 30 PSI models also available on request.



Samco Medical

Website: www.samcomedical.com,