



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Infant Weighing Scale - Beam Type

**Product Code :**  
SAMCOMEDICAL170002A020



**Description :**

Infant Weighing Scale - Beam Type Beam type (mechanical) scale for infants, 16 kg x 10 gm. Measuring range up to approx 16 kg. Minimum graduation, 10 g. Readout in kg only. Easy readable in low light working situations. Reading time max 15 seconds. Includes removable tray for infant. With easy zero adjustment possibility. Design allows rough handling. Smooth surface/finishing allows for easy cleaning/disinfection. All vital parts made of rust proof materials. Easy maintenance and repair in low tech settings. Splash proof and shock resistant light-weight body. Materials; body, metal (white baked enamel) or plastic Sturdy design allows heavy duty use.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),