



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Weighing Scale, Mechanical Column Type For Adults

**Product Code :**  
SAMCOMEDICAL170002A038



**Description :**

Weighing Scale, Mechanical Column Type For Adults Eye level sliding weights Height measuring rod, range 75 to 200cm / division 5mm Large non-slip platform, Robust manufacture, corrosion protected



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),