



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Weighing Scale, Standing, Adult, Digital with Height Rod

Product Code :

SAMCOMEDICAL170002A039

**Description :**

Weighing Scale, Standing, Adult, Digital with Height Rod Capacity: 300 Kg X 0.1Kg Accuracy: $\pm 150g$ Weight Unit: Kg/lb LCD Display: 1.0" LCD digital Display Dimension: Platform Size: 310 x 310 x 54mm, Column Height: 900mm Key Function: Hold/BMI, Unit/Weight, On/Zero/Off, Tare Power Supply: 1.5V 300mA adapter 1.5V AAA battery *6 Operation Temperature: 5°C / 35°C Height Rod: 60-205cm / 6ft 9inch



Samco Medical

Website: www.samcomedical.com,