



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Height Measuring Scale (Infant) / Infantometer

**Product Code :**  
SAMCOMEDICAL170002A040



**Description :**

Height Measuring Scale (Infant) / Infantometer Sleek broad base with one sliding side adjustable as per length of baby. Dual scale for direct reading in cms / inches upto 90 cm. Folding slides for easy storage. Baby can comfortably lay onto the chest of the scale



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),