



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Bathroom Scale electronics flat with tare parent model

Product Code :

SAMCOMEDICAL170002A045

**Description :**

Bathroom Scale electronics flat with tare parent model Maximum Capacity: 300 kg / 660 lbs. Graduation: 0.1 kg / 0.2 lb. Accuracy: $\pm 1/2$ 150g Resolution: 1:3000 Tare Range: Full Range Repeatability(STD dev): ± 1 d Linearity: ± 3 d Unit of Measure: Kg/lb Key Function: Hold/Unit/Weight, On/Zero/Off, Tare Stabilisation Time:1-2Sec Operating Temperature: $5\frac{1}{2}^{\circ}\text{C}$ - $35\frac{1}{2}^{\circ}\text{C}$ or $13\frac{1}{2}^{\circ}\text{F}$ - $95\frac{1}{2}^{\circ}\text{F}$ Power Supply: 9V battery or AC adaptor Calibration: Automatic External Calibration Mass: 300KG Indicator Display:1.0"LCD digital display Dimension(w x d x h): 15.5 x12.5x3.3 inch 39.5 x 32 x 8.5 cm Gross Weight : 7.7lb/3.5 Kg



Samco Medical

Website: www.samcomedical.com,