



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Vacuum Extractor Set

Product Code :
SAMCOMEDICAL180001A007



Description :

Vacuum Extractor Set Vacuum Extractor Set, Malmstorm type, manual operated. Complete with 3 S.S. Cups, (40, 50, 60mm). PVC pressure tube. Traction Handle. Suction Pump Glass jar in between top and bottom metal plates fitted with vacuum gauge and control valve. Packed in plastic case.



Samco Medical

Website: www.samcomedical.com,