



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Vacuum Extractor - Manual

**Product Code :**  
SAMCOMEDICAL180001A009



**Description :**

Vacuum Extractor - Manual Vacuum Extractor Set, Bird Type, manual operated. Complete with 3 stainless steel cups (40, 50, 60mm). Autoclavable Silicone pressure tube. Traction handle Suction pump, glass jar fitted between M.S. Basket with vacuum gauge and control valve



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),