



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Vacuum Extractor Electric with S.S. Cups

Product Code :
SAMCOMEDICAL180001A011



Description :

Vacuum Extractor Electric with S.S. Cups Housing: ABS Plastic Moulded Capacity: -575mmHg $\frac{1}{2}$ 10 , 45 Ltrs/min Pump Type: Diaphragm Pump Jars: 1 x 2 Ltrs Polycarbonate Jar Filter: Bacterial Filter Autoclavable/Reusable Tubing: 8mmID x 2mtr. (Silicone) Vacuum Gauge: 6.25 cm dia 0-760mmHg calibration. Power: 220 V AC, 50 / 60 Hz, 90Watt (110V on request) Noise Level: 55 dB A $\frac{1}{2}$ 3 Dimension & Weight: 43 x 22 x 37 cms, 8 Kg. Overflow Safety: Mechanical / Electronic Acoustic Alarm Supplied with 3 S.S. cups (40, 50, & 60mmdia.)



Samco Medical

Website: www.samcomedical.com,