



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Ideal Delivery Aid for Obstetricians : Safe for Mother / baby

Product Code :

SAMCOMEDICAL180001A012

**Description :**

Ideal Delivery Aid for Obstetricians : Safe for Mother / baby Housing: ABS Plastic Moulded Capacity: -600 mm Hg $\frac{1}{2}$ 10 , 12 Ltrs/ min Pump Type: Piston Pump Jars: 1 x 1 Ltr Polycarbonate Jar Filter: Bacterial Filter Autoclavable/Reusable Tubing: 8 mm ID x 2 mtr. (PVC) Vacuum Gauge: 5 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 / 60 Hz, 90Watt (110V on request) Noise Level: 55 dB A $\frac{1}{2}$ 3 Dimension & Weight: 30 x 14 x 33 cms, 4.5 Kg Overflow Safety: Mechanical / Electronic Acoustic Alarm Venteuse cup: Silicone - 50, 60 & 70mm : minimizes risk of trauma and cephalohematoma of fetal cranium



Samco Medical

Website: www.samcomedical.com,