



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hospital Plasma Thawing Bath

Product Code :
SAMCOMEDICAL230001A006



Description :

Hospital Plasma Thawing Bath - Plasma Thawing Bath, is a tabletop model, for quick thawing of plasma at 37 oC. In this, the temperature is controlled by micro controller based temperature controller. Our range is designed to safely and reliably thaw fresh frozen plasma (FFP) for the recovery of cryoprecipitated antihemophilic factor (AHF) and in this, the temperature is preset at 4 oC and maintained with 0.2 oC stability Hospital Plasma Thawing Bath, Blood Bank Equipments Exporters, Blood Bank Equipments Manufacturers



Samco Medical

Website: www.samcomedical.com,