



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Aluminium Case for Midwifery Kit

Product Code :
SAMCOMEDICAL180001A026



Description :

Aluminium Case for Midwifery Kit Case, Aluminium, for Midwife Kit , approx 305 x 216 x150 mm (12" x 8-1/2" x 6") with removable tight fitting lid collapsible carrying handle on cover, and one sturdy spring loaded clamp at each small end as a locking device All corners fully rounded and all edges smooth finished.



Samco Medical

Website: www.samcomedical.com,