



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Hospital Blood Bag Tube Sealer

**Product Code :**  
SAMCOMEDICAL230001A009



**Description :**

Hospital Blood Bag Tube Sealer, which is a compact equipment to seal the blood bag pilot tube by ratio frequency sealing system. The sealing can be achieved in less than 1.5 seconds per tube by simple operation. Apart from this, the sealed tube can be easily separated by pulling the sealing section by both sides. Hospital Blood Bag Tube Sealer, Blood Bank Equipments Exporters, Blood Bank Equipments Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),