



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Suction Autoclave Jar

Product Code :
SAMCOMEDICAL150008A001

Description :

PC Suction Jar with Lid Capacity of Jar 2Lt.: Autoclave able Safety Jar with ABS lid. The jars graduated at 50 ml intervals. The jars are capable for repeated sterilization at 120 oC. Capacity of jars are 2.0 liters. The Jars are made of Polycarbonate without any pronounced tint and any visual defects. The lid and its connections provide a leak proof assembly. Jar & lid are comes with overflow safety arrangement. Rubber Gaskets satisfy the aging test at 70 oC Suction Autoclave Jar, Parts Suction Units Exporters, Parts Suction Units Manufacturers



Samco Medical

Website: www.samcomedical.com,