



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
High Vacuum Suction Machine

**Product Code :**  
SAMCOMEDICAL150003A001

**Description :**

Non Corrosive Metal Cabinet with Stainless Steel Top. Rotary Vane Type Vacuum Pump with -730 (+-10 mm Hg. Free air displacement 40 - 45 LPM, Noise level of suction unit is less than 50 dBA, Electrical Requirement 220-230VAC, 50Hz, 1Phase. Heavy duty HN-65 Castors. 0.50 H.P, 380 Watt, 1440 RPM Electric Motor. 100 mm Calibrated Vacuum Gauge. Non collapsible Suction tubing. 2 X 2.0 Ltrs. Polycarbonate jars with overflow safety. Bacterial filter fitted on top High Vacuum Suction Machine, High Vacuum Suction Exporters, High Vacuum Suction Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),