



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Led Shadowless Operation Lamp

**Product Code :**  
SAMCOMEDICAL200001A050



**Description :**

Led Shadowless Operation Lamp Technical Specifications: Illumination : 160000 LUX Color temperature : 4300i¿½500K Color indication index (Ra) : 95 Power of LED bead : N128x1W Spot diameter d50 : 130mm-350mm Illumination depth : 140cm Illumination-adjustment range : 10%~100% Ultraviolet irradiance with wavelength under 400nm : 0.002W/mi¿½ Focusing control : Manual Perfect for endoscopy surgery Power supply : 220V, 50Hz /110V, 60Hz LED average service life : 60000hrs Lowest height of installation : 2800mm Available in ceiling, floor and wall configurations



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),