



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ceiling Shadowless Operation Lamp Single Reflector

Product Code :
SAMCOMEDICAL200001A053



Description :

Ceiling Shadowless Operation Lamp Single Reflector Illumination Intensity : 80000-150000 Lux Color Temperature : 3500-5000 K Spot Diameter : 100-250 mm Illumination Depth : 600 mm Height Adjustment : 800 mm Lamphead Diameter : 700 mm Brightness Control : 9-level dimmable Temperature Rise of Field : ≤ 2 ($1\frac{1}{2}^{\circ}\text{C}$) Power Supply : 220V / 50Hz Power Input : 200 VA Specification of Bulb : 150 Watt / 24 Volt Minimum Height for Installation : 2700 mm Automatic Switchover Between the Key and Backup Bulbs : Yes



Samco Medical

Website: www.samcomedical.com,