



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ceiling Shadowless OT Light Single Dome (6 Bulbs)

Product Code :
SAMCOMEDICAL200001A056



Description :

Ceiling Shadowless OT Light Single Dome (6 Bulbs) Illumination Intensity : 100000~160000 Lux Color Temperature : 3500-5000 K Spot Diameter : 120-250mm Height Adjustment : 700mm Lamphead Diameter : 800mm Brightness Control : 9-level dimmable Temperature Rise of Field : ± 2 ($\pm 1/2^\circ\text{C}$) Power Supply : 220V / 50Hz Power Input : 400 VA Specification of Bulb : 50Watt / 24 Volt Minimum Height for Installation : 2750mm



Samco Medical

Website: www.samcomedical.com,