



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
High Vacuum Suction LIPO

**Product Code :**  
SAMCOMEDICAL150003A003



**Description :**

Complete S.S Body to prevent from rust. Heavy Duty 1 hp motorized pump for fast vacuum built  $\frac{1}{2}$  up. Absolutely Noiseless. Self lubricating system for trouble free service. Finely regulable 0-750mm Hg vacuum, Easy to read off. Over Flow Safety Device: over  $\frac{1}{2}$  flow check valve and removable filter to protect the pump from moisture. Non collapsible tubing Max Vacuum: 730  $\frac{1}{2}$  10mm Hg regulable through Vacuum Control Knob. Suction Capacity: 40 to 45Lt per Min Vacuum Gauge Bourden Type 4 0-750mm Hg calibration Suction Bottle: Glass jars or 2 x 2000cc Poly carbonate jars ( UNBREAK ABLE ) with over flow check valve (float type) Motor: 1 hp motor Vacuum Pump: Rotary type ( double pump ) . Power Source: AC 220V 50 Hz. Single Phase. Overall Size: L 380mm x W 390mm x H 825mm (Approx.) Finish: Powder coated. Weight: 34 kg ( approx ) High Vacuum Suction LIPO, High Vacuum Suction Exporters, High Vacuum Suction Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),