



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Ceiling Shadowless Operation Lamp Single Dome (7 Bulbs)

Product Code :

SAMCOMEDICAL200001A057

**Description :**

Ceiling Shadowless Operation Lamp Single Dome (7 Bulbs) Illumination Intensity : 42000-70000 Lux Color Temperature : 3500-5000 K Spot Diameter : 100-200mm Height Adjustment : 440mm Lamphead Diameter : 700mm Power Supply : 220V / 50Hz Power Input : 210 VA



Samco Medical

Website: www.samcomedical.com,