



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ceiling Shadowless Operation Lamp Single Dome (9 Bulbs)

Product Code :
SAMCOMEDICAL200001A058



Description :

Ceiling Shadowless Operation Lamp Single Dome (9 Bulbs) Illumination Intensity : 54000-90000 Lux Color Temperature : 3500-5000 K Spot Diameter : 100-200mm Height Adjustment : 440mm Lamphead Diameter : 900mm Power Supply : 220V / 50Hz Power Input : 270 VA



Samco Medical

Website: www.samcomedical.com,