



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Ceiling Mounted Shadowless Operation Lamp With Single Reflector (With Focus Control)

**Product Code :**  
SAMCOMEDICAL200001A060



**Description :**

Ceiling Mounted Shadowless Operation Lamp With Single Reflector (With Focus Control) Illumination Intensity : 30000-60000 Lux Color Temperature : 3000-4500 K Spot Diameter : 80-150mm Height Adjustment : 700mm Lamphead Diameter : 195mm Power Supply : 200V / 50Hz Power Input : 75 VA Specification of Halogen Bulb : 50Watt / 24 Volt Brightness Control : Digital Control



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),