



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ceiling Mounted Led Surgical Light Dual Arm - 4+3 Lamp

Product Code :
SAMCOMEDICAL200001A062



Description :

Ceiling Mounted Led Surgical Light Dual Arm - 4+3 Lamp Diameter of lamp head : Main lamp, -700mm, Subsidiary, -600mm Illumination : Main lamp, 80000 lx ~ 160000 lx : Subsidiary, 70000 lx ~ 160000 lx Color temperature : 3000k ~ 6000k Color rendering index : 92 ~ 100 Diameter of Spot : D10 Main lamp 140mm, Subsidiary 120mm : D50 Main lamp 80mm, Subsidiary 70mm Depth of illumination : = 600mm Adjustment of Height : 800mm Brightness control : 9 steps Total radiation illumination : = 500W/M2 Power supply : A.C100 ~ 220V 50 HZ Power consumption : 120VA Power of bulb : 3 WLED Number of LED Bulbs : Main lamp, 28 : Subsidiary, 21 Minimum installation height : 3m



Samco Medical

Website: www.samcomedical.com,