



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Suction Tube

**Product Code :**

SAMCOMEDICAL150008A002

**Description :**

Suction Tube: Silicon or PVC transparent tubing with soft sealing grips. Retain elasticity even after prolonged use. ID =10mm & wall thickness = 2.5 mm. Excellent flexible & heat resistant. Odor less & repeatedly autoclaveable. Transparency of tube allows visibility of flowing fluids Suction Tube, Parts Suction Units Exporters, Parts Suction Units Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),