



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ceiling Shadowless Operation Lamp Dual Arm (6+3 Bulbs)

Product Code :
SAMCOMEDICAL200001A067



Description :

Ceiling Shadowless Operation Lamp Dual Arm (6+3 Bulbs) Illumination Intensity : 100000~160000
60000~100000 Lux Color Temperature : 3500-5000 K Spot Diameter : 120-250mm Height Adjustment : 700 /
800mm Lamphead Diameter : 800 / 550mm Brightness Control : 9-level dimmable Temperature Rise of Field :
?2 (1½C) Power Supply : 220V / 50Hz Power Input : 600 VA Specification of Bulb : 50Watt / 24 Volt Minimum
Height for Installation : 2850mm



Samco Medical

Website: www.samcomedical.com,