



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ceiling Shadowless Operation Lamp Dual Arm (5+12
Bulbs)

Product Code :
SAMCOMEDICAL200001A068



Description :

Ceiling Shadowless Operation Lamp Dual Arm (5+12 Bulbs) Illumination Intensity : 50000-120000 Lux Color Temperature : 3500-5000 K Spot Diameter : 100-200mm Height Adjustment : 260-440mm Lamphead Diameter : 500 / 900mm Power Supply : 220V / 50Hz Power Input : 510 VA



Samco Medical

Website: www.samcomedical.com,