



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Shadowless Operation Lamp With Single Reflector Stand
Model With Castors (With Focus Control)

Product Code :

SAMCOMEDICAL200001A071

**Description :**

Shadowless Operation Lamp With Single Reflector Stand Model With Castors (With Focus Control) Illumination Intensity : 30000-60000 Lux Color Temperature : 3000-4500 K Spot Diameter : 80-150mm Height Adjustment : 700mm Lamphead Diameter : 195mm Power Supply : 200V / 50Hz Power Input : 75 VA Specification of Halogen Bulb : 50Watt / 24 Volt Brightness Control : Digital Control



Samco Medical

Website: www.samcomedical.com,